



M1



M2



M3



M4



M5



M6



M7



M8

## Mini 系列可选模块

### 可更换模块

型号	试管直径	试管数	型号	试管直径	试管数
M1	0.2ml 离心管	32	M5	比色皿	8
M2	0.5ml 离心管	24	M6	15ml 离心管 (无法合盖)	6
M3	1.5ml 离心管	15	M7	50ml 离心管 (无法合盖)	2
M4	2.0ml 离心管	15	M8	5ml 离心管 (无法合盖)	12